



Colloidal Oatmeal System

Nitrile Exam Gloves Powder Free, Standard Cuff

COATS® (an acronym for colloidal oatmeal system) is a patented and unique nitrile glove technology. COATS® utilises the powerful benefits of all-natural oats, an FDA-recognised skin protectant, as a coating that forms a natural, moisturising barrier between the glove and skin. This acts as a preventative measure against skin irritation, and eliminates many of the uncomfortable and irritating conditions experienced when wearing normal gloves. Users who suffer from dry and itchy skin due to constant hand washing and glove usage can now rely on COATS® to soothe and nurture the skin, and protect their hands while they work



Colloidal Oatmeal Content

	COATS® Nitrile			
Length (mm)				
	<u>≥</u> 230			
Thickness Measurements (mm)				
Palm (centre of Palm)	0.07 ± 0.02			
Finger (13mm \pm 3mm from tip)	0.09 ± 0.02			
Physical Properties	Before Ageing	After Ageing		
Tensile Strength (MPa)	≥ 18	≥ 16		
Elongation (%)	≥ 500	≥ 400		
Inspection Levels & AQL	Inspection Level	AQL		
Watertightness	G1	1.5		
Physical Dimensions	S2	4.0		
Physical Properties	S2	4.0		
Visual Inspection (Major)	S4	2.5		
Visual Inspection (Minor)	S4	4.0		
Particulate Residue	N = 5	≤ 2mg/glove		

Chemotherapy Drugs and Concentration (Tested for Resistance to Permeation by Chemotherapy Drugs as per ASTM D6978-05)	Minimum Breakthrough Detection Time (minutes)
Carmustine (BCNU), 3.3mg/ml (3,300 ppm)	Not recommended
Cisplatin, 1.0mg/ml (1,000 ppm)	>240 minutes
Cyclophosphamide (Cytoxan), 20.0mg/ml (20,000 ppm)	>240 minutes
Dacarbazine (DTIC), 10.0mg/ml (10,000 ppm)	>240 minutes
Doxorubicin Hydrochloride, 2.0mg/ml (2,000 ppm)	>240 minutes
Etoposide (Toposar), 20.00mg/ml (20,000 ppm)	>240 minutes
Fluorouracil, 50.0mg/ml (50,000 ppm)	>240 minutes
Methotrexate, 25.0mg/ml (25,000 ppm)	>240 minutes
Mitomycin C, 0.5mg/ml (500 ppm)	>240 minutes
Paclitaxel (Taxol), 6.0mg/ml (6,000 ppm)	>240 minutes
Thiotepa, 10.0mg/ml (10,000 ppm)	Not recommended
Vincristine Sulfate, 1.0mg/ml (1,000 ppm)	>240 minutes

WARNING: Gloves used for protection against chemotherapy drug exposure should be selected specifically for the type of chemicals being used. Due to the variety and concentration of chemotherapy drugs used in treatments, the resistance table shown does neither warrant nor imply the safe use of the gloves against chemotherapy drugs resistance in every case. The safe use of gloves in chemotherapytreatment is solely the decision of clinicians authorised to make such decision.

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Standard cuff
- Dawn blue colour

PACKAGING

100 gloves per box (XS-L) 90 gloves per box (XL) 10 boxes per carton

REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k), MDD 93/42/EEC, REACH, EC 10/2011, EC 1935/2004

STANDARDS

ASTM D6319, ASTM D412, ASTM D573, ASTM D5151, ASTM D6124, EN 455 part 1, 2, 3 & 4, EN 1186, EN 13130, CEN/TS 14234

PATENTS

Patent 7,691,436; Patent 7,718,240; Patent 7,740,622; Patent 8,075,965; Patent 8,458,818

MANUFACTURING ACCREDITATIONS

ISO 9001 ISO 13485 EN ISO 13485 ISO 14001 OHSAS 18001



COATS® Colloidal Oatmeal Coated Nitrile Powder Free 2.5 Mil

ASTM D3578

Physical Dimensions		
Glove Length (mm)	≥ 230	
Palm Thickness (mm)	0.07 ± 0.02	
Finger Thickness (mm)	0.09 ± 0.02	
Physical Properties		
Test	Before Aging After Agi	
Tensile strength (MPa)	≥ 18.0	≥ 16.0
Elongation (%)	≥ 500	≥ 400

EN 455

Physical Dimensions			
Median glove length (mm) ≥ 240		40	
Median palm thickness (mm) 0.07 ± 0.02		0.02	
Median finger thickness (mm)	0.09 ± 0.02		
Physical Properties			
Test	Before Aging	After Aging	
Median Force at break (N)	≥ 6 ≥ 6		



Regulatory Compliance

FDA 510(k), MDD 93/42/EEC, REACH, ROHS Directive 2002/95/EC, EC 10/2011, EC 1935/2004, PPE 89/686/EEC

Standards

ASTM D6319, ASTM 6978, EN455 part 1, 2, 3 & 4, EN 1186, EN 13130, CEN/TS 14234, EN 420, EN 374 part 1, 2 & 3

Classification

Class I (FDA), Class I (MDD 93/42/EEC), Category 3 (BfR XXI), Category III (PPE 89/686/EEC)

Patent

7,691,436; 7,718,240; 7,740,622; 8,075,965; 8,458,818


Application Settings

Low risk - medical, dental, procedures, chemotherapy drugs, pathology lab and food handling. Coated with FDA recognised skin protectant. Clinically proven to help protect and moisturise your skin from dry and irritated skin from prolonged glove use and hand wash.

Colour

Dawn blue, white









COMMON NAME (USED ON LABEL) Nitrile Powder Free Examination Gloves	CHEMICAL FAMILY Carboxylated Butadiene Acrylonitrile Polymer Latex
APPLICATION	TRADENAME & SYNONYM
Medical and Dental	GLOVEON COATS NITRILE (CTS38)
	NITRILE POWDER FREE EXAMINATION GLOVES COATS

SECTION 2: HAZARDOUS INGREDIENTS					
HAZARDOUS COMPONENT	CAS #	%(WT)	TLV	PEL	
N/A	N/A	N/A	N/A	N/A	

PEL: Permissible Exposure Limit established by Occupational Safety and Health Administration (OSHA).

TLV: Threshold Limit Value established by the American Conference of Governmental Industrial Hygienists, 1987-1988.

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

CHEMICAL COMPOSITION

All chemicals used are non-toxic/ non-hazardous.

Butadiene-Acrylonitrile Latex, Sodium Dodecylbenzenesulfonate, Sulphur, Zinc Oxide, Zinc Di-nbutyldithiocarbamate, Titanium Dioxide, Paraffin Wax Emulsion

Coating Ingredient

Colloidal Oatmeal & Constituents, Sodium Benzoate, Processing Aid

SECTION 4: FIRST AID MEASURE

If reaction in the form of skin irritation is noticed, remove gloves immediately and wash affected part with saline water. If there is no relief, seek medical reactions.

SECTION 5: FIRE FIGHTING MEASURE FLASHPOINT AUTOIGNITION TEMPERATURE FLAMMABLE LIMITS IN AIR

EXTINGUISHING MEDIA

Chemical foam and dry chemical may be used.

FIRE-FIGHTING PROCEDURES

Use standard procedures for combustion material fires, including approved self-contained breathing apparatus.

FIRE AND EXPLOSION HAZARDS

No fire or explosion hazards are associated with these products. They will melt at elevated temperatures.

SECTION 6: ACCIDENTAL RELEASE MEASURES

BIOCOMPATABILITY

The chemical formulation of the gloves and surface lubrication materials does not contain any substances normally known to be harmful to the user or to any person with whom the gloves come into contact.

MEDICAL CONDITIONS GENERALLY AGGRAVATED BY EXPOSURE

Nitrile Powder Free Gloves are not expected to cause any adverse health effects.

SECTION 7: HANDLING AND STORAGE

PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE

Store in a dry, cool and ventilated area. Do not store above 104 °F (40 °C). Shield open box from direct sunlight, fluorescent lighting and x-rays. Improper storage will decrease usable life.



SECTION 8: EXPOSURE CONTROLS/ PERSONAL PROTECTION **EYE PROTECTION SKIN PROTECTION** Not necessary under conditions of intended use. Not necessary under conditions of intended use. RESPIRATORY PROTECTION VENTILATION Not necessary under conditions of intended use. Not necessary under conditions of intended use.

STEPS TO BE TAKEN IN CASE MATERIAL IS LEAKED OR SPILLED

These products are solid articles and are not subject to leaks or spills.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE/ ODOR

Ambidextrous, Beaded Cuff, Micro-textured, Chlorinated, Powder Free, Coated with Colloidal Oatmeal USP Skin Protectant, Dawn Blue.

DIMENSION	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
Length (mm)	S	Mi	nimum 230 (same fo	or all)	ž.
Width (mm)	76 ± 4	86 ± 4	98 ± 4	107 ± 4	115 ± 4
THICKNESS (mm	- SINGLE WALL	MEASUREMENT (same for all)	-	•
Control of the state of the sta			Venney Control		

Finger (mm) 0.09 ± 0.02 Palm (mm) 0.07 ± 0.02

TENSILE PROPERTIES	UNAGED	AGED	
Tensile Strength (Mpa)	Min. 18.0 MPa	Min. 16.0 MPa	
Ultimate Elongation (%)	Min. 500%	Min. 400%	

SECTION 10: STABILITY AND REACTIVITY

BOILING POINT	VAPOR PRESSURE (mm Hg)	VAPOR DENSITY (air=1)
N/A	N/A	N/A
SPECIFIC GRAVITY (water=1) N/A	SOLUBILITY IN WATER Insoluble	% VOLATILE BY VOLUME N/A

EVAPORATION RATE	VISCOSITY
N/A	N/A

SECTION 11: TOXICOLOGICAL INFORMATION

STABILITY **CONDITIONS TO AVOID** Stable. Does not apply.

INCOMPATABILITY (MATERIALS TO AVOID)

High polar solvent like methyl ethyl ketone, acetone.

HAZARDOUS DECOMPOSITION PRODUCTS

In a fire, these products may produce a black smoke. Carbon Dioxide, Carbon Monoxide, Oxides of Nitrogen, aromatic/aliphatic hydrocarbons.

HAZARDOUS POLYMERIZATION

Will not occur.

SECTION 12: ECOLOGICAL INFORMATION

N/A

SECTION 13: DISPOSAL CONSIDERATION

WASTE DISPOSAL METHOD

Consult current local, state and federal regulations for proper disposal methods.

SECTION 14: TRANSPORT INFORMATION

N/A

SECTION 15: REGULATORY INFORMATION

N/A

SECTION 16: OTHER INFORMATION

RECOMMENDED USE AND RESTRICTION

The Nitrile Powder Free Gloves is a Single Use device.



The Brand

Leadership | Certifications | Global Locations

Certifications

Certifications

Gloveon's quality standards, management systems and exemplary regulatory compliance, all contribute to the global success of the company. Our capabilities have been assessed and certified by the following international governing bodies.



Management Service ISO 9001:2015



America ISO 13485:2016



EN ISO 13485:2016



Confirmation Letter for GMP Audit



EC Certificate



ISO 14001:2015



UL Certification



ISEGA Food Contact Test Certification (German)



Registration Certificate for Medical Device



NFPA Certification



510(k) Approval



PPE Cert















Report No. : CRSSA/02645/18

TEST REPORT

Product Description Country of Origin Size Quantity Tested Test Conducted Test Method Testing Period Powder Free Nitrile Examination Gloves Malaysia Medium 200 pieces Freedom from holes EN455 Part 1:2000 02 Mar 2018 – 08 Mar 2018

Based on submitted samples the following results obtain

: Within AQL

Note: Upon Customer's request, this report has been issued in more than one original. Only the first original is a legally binding document and may be used for any legal purpose, including payment. (Original 1-3)

SGS (MALAYSIA) SDN. BHD.

Chee. CHEE TUCK CHOON B.Sc. MMIC SECTION HEAD

zer restatient for serves days (perhaltees noce samples or more morms cary.

SGS (Malahysis) SAG: BMd.

(Company No. 10871-T)

4+6 (03) 5481 8282 f +6 (03) 5481 8215 www.sgs.com





Hartalega Sdn. Bhd. Nurul Kong Quality Assurance Senior Manager No. 7, Kawasan Perusahaan Suria Bestari Jaya, 45600 My

Re: K180505
Trade/Device Name: Nitrile Powder Free Examination Glove with Colloidal Oatmeal USP with LowDermatitis Potential Claim and Tested For Use with Chemotherapy Drugs (White)


Regulation Number: 21 CFR 880.6250 Regulation Number: 21 CFR 880.6250 Regulation Name: Patient Examination Glove Regulatory Class: Class I Product Code: LZA, LZC Dated: May 15, 2018 Received: May 17, 2018

Dear Nurul Kong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced We have reviewed your Section \$10(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act michael requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations. Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to registration and listing [21 CFR part 807); labeling [21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good



Page 2 - Nurul Kong

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/MedicalDevices/Safety/BeportaPublicaldefault.htm.

For comprehensive regulatory information about medical dockies and realizion-emitting products, including information about labeling regulations, please see Device Advice (https://www.fds.nov.MedicalDevices/Device/Devic

Clarence W. Murray

Enclosure

TMENT OF HEALTHAND HUMAN SE Indications for Use

K180505

r Name
Pornder Free Examination Glove with Colloidal Outmail USP with Low Demutitis Potential Coim and Tested for Use with shoupy Brups (White)

Initiation for Use (Charoba)
Nitrik Powder Free Exemination Glove with Calloidd Outmed USP with Low Demantis Potential Claim and Tosted for De with Chemishenpy Drugs (White) is a non-stelle deposable derive intended for medical purpose that is wom on the construct least to prevent contensionion between patient and constinct. It is also tend to be used against Chemisthenpy

The list of Chemothempy Drugs tested (with breakthrough times) is as belo


17.2 minutes >240 minutes Carmanine (3.3 mg/ml)
Caplain (1.6mg/ml)
Cyclopharine (2.0 mg/ml)
Cyclopharphanide (2.0 mg/ml)
Dascarbaine (0.0 mg/ml)
Dascarbisine Hydrochloride (2.0 mg/ml)
Daponide (2.0 mg/ml)
Flantouracii (30 mg/ml)

CONTINUE ON A SEPARATE PAGE F NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1925.
"DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW:

Page 1of 1






Hartalega





The voluntary certification mark with the statement "Type tested" is issued for products and components. The certification mark demonstrates that the

Familiar from hor tools and toys, the indicates that a pr safety according to





LOCATIONS CONTACT US Q

















