

FJT-200 Series Flocculation Jar Testers

- Digital display clearly indicates speed of paddles
- Save time with two preset programs
- BioCote® antimicrobial protection is built in for a cleaner, more hygienic product
- Built-in illumination and dark background for easy viewing of samples
- Two or six stirring positions are ideal for small-scale testing
- 3-year warranty



FJT-200D Flocculator
(beakers not included)



FJT-200D-6 Flocculator
(beakers not included)

Specifications

Model no.	FJT-200D	FJT-200D-6
No. of samples	2	6
Speed range	25 to 250 rpm	25 to 250 rpm
Timer	No	0 to 99 minutes
Pre-set programs	2	2
Controls keypad	Dial	Touch sensitive
Digital displays	LCD	LED
Dimensions (W x H x D)	32 x 50 x 21 cm (12.6 x 19.7 x 8.3 in)	75 x 46 x 21 cm (29.5 x 18.1 x 8.3 in)
Electricity supply	230 VAC, 50/60 Hz, 100 W	230 VAC, 50 Hz, 200 W
IP rating	IP31	IP31
Warranty	3 years	3 years

FJT-200 Series Flocculation Jar Testers

Our FJT-200 flocculation jar testers are designed to test the efficiency of flocculation or precipitation reagents on a small pilot-scale. The single motor drives all positions at a constant stirring speed. Digital display clearly indicates the speed of the rotational stainless steel paddles.

Two preset programs allow routine speed and speed/time parameters to be input and then selected by simply pressing one button. Parameters are stored in the memory and can be recalled or altered at any time.

Sample monitoring is assisted by diffused base illumination and black background. Both testers accept 1-L beakers (not included).

FJT-200D Flocculator accepts two samples (without timer function). Tester features speed control knob and LCD screen. The case incorporates a handle for easy transportation.

FJT-200D-6 Flocculator accommodates up to six samples and is specially designed for repeatable conditions between samples and from run to run (timer function 0 to 99 minutes). Tester features a completely flat keypad that is easy to clean, and touch-sensitive controls and LED displays for time and speed. Order optional coagulant doser to dispense coagulant chemicals simultaneously.

Ordering Information

Description	Item no.	Series no.	Model no.	Legacy SKU
Flocculation Jar Tester, 2 positions	99551-00	FJT-200	FJT-200D	SW5
Flocculation Jar Tester, 6 positions	99551-05	FJT-200	FJT-200D-6	SW6

Accessories

Description	Item no.	Series no.	Model no.	Legacy SKU
Simultaneous coagulant doser for 6-position model only	99551-99	FJT-200	N/A	SW6/2

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1 800 363 5900
China: +86 21 5109 9909
France: +33 (0)1 48 63 78 00
Germany: +49 (0) 9377 92030

India: +1 800 266 1244
Italy: +39 (0)2 84349215
UK: +44 (0) 1480 272279
All other countries: +1 847 549 7600

Cole-Parmer®